

AD-A038 836

ARMY MATERIALS AND MECHANICS RESEARCH CENTER WATERTO--ETC F/G 14/5
A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION--ETC(U)
JAN 77 B H STRAUSS

UNCLASSIFIED

AMMRC-TR-77-1

NL

| OF |
AD
A038836



END

DATE
FILMED
5-77

AD A 038836

AMMRC TR 77-1


12

AD

A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC ANALYSIS

BERNARD H. STRAUSS
POLYMER AND CHEMISTRY DIVISION

January 1977



Approved for public release; distribution unlimited.

AU NO. _____
DDC FILE COPY

ARMY MATERIALS AND MECHANICS RESEARCH CENTER
Watertown, Massachusetts 02172

0A

AMMRC TR 77-1

A COMPUTER METHOD
FOR THE CALIBRATION OF
A PHOTOGRAPHIC EMULSION
FOR SPECTROGRAPHIC ANALYSIS

The findings in this report are not to be construed as an official Department of the Army position, unless so designated by other authorized documents.

Mention of any trade names or manufacturers in this report shall not be construed as advertising nor as an official indorsement or approval of such products or companies by the United States Government.

DISPOSITION INSTRUCTIONS

Destroy this report when it is no longer needed.
Do not return it to the originator.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER AMMRC-TR-77-1	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC ANALYSIS	5. TYPE OF REPORT & PERIOD COVERED Final Report	
7. AUTHOR(s) Bernard H. Strauss	8. CONTRACT OR GRANT NUMBER(s)	
9. PERFORMING ORGANIZATION NAME AND ADDRESS Army Materials and Mechanics Research Center Watertown, Massachusetts 02172 DRXMR-RA	10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS D/A Project: A1-6-P6350-AW-AW AMCMS Code: 5397.OM.6350	
11. CONTROLLING OFFICE NAME AND ADDRESS U. S. Army Materiel Development and Readiness Command, Alexandria, Virginia 22333	12. REPORT DATE January 1977	13. NUMBER OF PAGES 20
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office) 1223p.	15. SECURITY CLASS. (of this report) Unclassified	
15a. DECLASSIFICATION/DOWNGRADING SCHEDULE		
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES This project has been accomplished as part of the U. S. Army Materials Testing Technology Program, which has for its objective the timely establishment of testing techniques, procedures or prototype equipment (in mechanical, chemical, or nondestructive testing) to insure efficient in- spection methods for materiel/material procured or maintained by DARCOM.		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Computer programs Glass Calibrating Emulsions Spectrographic analysis		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) (SEE REVERSE SIDE)		

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

403 105

LB

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

Block No. 20

ABSTRACT

↘
A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.
↑

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

CONTENTS

	Page
INTRODUCTION.	1
PROCEDURE FOR EMULSION PLATE CALIBRATION.	2
HEWLETT-PACKARD CALCULATOR METHOD	4
UNIVAC 1108 COMPUTER METHOD	7
CONCLUSION.	8
ACKNOWLEDGMENTS	8
APPENDIX A. SPECTROGRAPHIC PARAMETERS.	9
APPENDIX B. HEWLETT-PACKARD CALCULATOR PROGRAM FOR GLASS PLATE EMULSION CALIBRATION.	10
APPENDIX C. UNIVAC 1108 PROGRAM.	15

NTIS	White Section	<input checked="" type="checkbox"/>
DDC	Buff Section	<input type="checkbox"/>
UNANNOUNCED		<input type="checkbox"/>
JUSTIFICATION		
BY _____		
DISTRIBUTION, AVAILABILITY CODES		
Dist.	AVAIL. AREA	SPECIAL
A		

INTRODUCTION

In the past ten years, many reports have been published concerning the application of computer programming for the calibration of photographic emulsions. Calculating boards¹ and rules² have essentially been replaced by readouts from large and small computers. Authors have generally chosen to use programs that follow graphical methods of calibration,³⁻⁸ although some, like Margoshes and Rasberry,⁹ take a different approach by assuming a linear relation between functions of the relative exposure and of the microphotometer readings. Most graphic methods studied appeared to be unduly complex; for example, data was generally converted to Seidel or "modified" Seidel transformations, the use of the principles involved in Kaiser transformations were employed, or dedicated computers were used.

Our laboratory has striven for a program that would follow our present, slightly modified ASTM method of calibration¹⁰ which is simple in principle. We sought a computer method that would plot preliminary and emulsion calibration curves so that a convenient check could be made on the original data at any time desired. By retaining the ability to have curves printed by the computer one can decide, for instance, whether more points are needed to establish the preliminary curve, as it must be remembered that the accuracy of the final readout of any computer program is limited by the accuracy, number, and selectiveness of the original dark/light data read from the microphotometer.

Using our program, dark/light data are inputted to obtain the emulsion calibration equation. This equation is then stored either in the memory of the calculator or on recording tape. To find a relative intensity, percent transmission readings (one from the element in question and the other from an internal standard) are entered in the computer. Troubleshooting the emulsion calibration is simplified by printing the preliminary curve with the original data superimposed on it along with the emulsion calibration curve.

Two similar programs were originated to follow our existing "hand emulsion calibration" method. One program uses time-sharing computer facilities and has the capability of printing charts or tables of percent transmission and corresponding relative intensities. The other program uses a Hewlett-Packard (H.P.) 9830-A calculator with an 8K memory storage. An X-Y plotter is attached to plot curves originated by the data from the computer.

1. KERBYSON, J. D. *A Calculating Board Incorporating a Linear Density Transform*. *Spectrochimica Acta*, v. 19, 1963, p. 1335-1341.
2. CHURCHILL, J. R. *Techniques of Quantitative Spectrographic Analysis*. *Ind. Eng. Chem. Anal. Ed.*, v. 16, 1944, p. 653-670.
3. HOLCOMBE, J. A., BRINKMAN, D. W., and SACHS, R. D. *Fortran-Based Photographic Emulsion Calibration Procedure for Use in Quantitative Spectrometry*. *Analytical Chemistry*, v. 47, March 1975, p. 441-447.
4. CARNEVALE, A., and LINCOLN, A. J. *A Computer Program to Determine Photographic Emulsion Calibration Curves*. *Spectrochimica Acta*, v. 24B, 1969, p. 313-323.
5. ANDERSON, J. W., and LINCOLN, A. J. *General Mathematical Approach to Emulsion Calibration in Optical Emission Spectroscopy*. *Applied Spectroscopy*, v. 22, 1968, p. 753-757.
6. JOYCE, C. S. *The Use of a Digital Computer for Photographic Photometry and Spectrochemical Computations*. *Canadian Spectroscopy*, v. 10, 1965, p. 33-37, 43.
7. GORDON, A., and GALLAGHER, A. K. *Automatic Data Processing for Photographic Photometry in Spectrographic Analysis*. Lewis Research Center, Cleveland Ohio, NASA TM X-1220, April 1966.
8. SHAW, D. M. *Principles of Machine Computation of Spectrochemical Analyses*. *Canadian Spectroscopy*, v. 10, 1965, p. 3-6.
9. MARGOSHES, M., and RASBERRY, S. *Application of Digital Computers in Spectrochemical Analysis—Computational Methods in Photographic Microphotometry*. *Spectrochimica Acta*, v. 24B, 1969, p. 497-513.
10. *Methods for Emission Spectrochemical Analysis*. American Society for Testing and Materials, 1971.

PROCEDURE FOR EMULSION PLATE CALIBRATION

Cassette, rotation, and position settings on the spectrograph were adjusted so that the wavelength region 2750 to 3450 Å could be photographed. A solution containing 50 mg/ml iron was sparked for a twenty- and forty-second exposure. Details concerning the parameters used for excitation and conditions set on the spectrograph are found in Appendix A. A seven-step filter, with segments in the ratio of 1:2, was placed between the spark and the entrance slit of the spectrograph to stepwise attenuate the exposure.

Using a densitometer, lines on the photographic plate in the wavelength region 3000 to 3200 Å were examined. Transmittances between 1.5% and 92% in the 66% and 33% portion of a line (the second and third step of the filter) were listed in a table labeled DARK percent T and LIGHT percent T. The preliminary curve (PC) is generated from data which consists of about 50 line pairs plotted and connected on linear graph paper. Paper 15" x 10" in size was labeled DARK percent T on the ordinate and LIGHT percent T on the abscissa. Formation of this curve is the first stage in the ASTM method.

The second stage of the method is the construction of the emulsion calibration or H and D curve and it is the equation of this curve that is ultimately used in relating percent transmission to relative intensity. To construct this curve a new table was established using the two parameters, percent T and NUMBER. A DARK percent T point was chosen on the preliminary curve as low in value as possible but above the intercept of the curve with the ordinate. A corresponding LIGHT percent T point was obtained and this value was written under percent T and 7 under NUMBER as in the table. Allowing the value of the LIGHT percent T to equal the numerical value of a DARK percent T point, a new LIGHT percent T value was obtained from the PC. In the table, 6 was placed under NUMBER and the new LIGHT percent T value recorded. This method was continued until LIGHT percent T values were obtained for the corresponding numbers 5, 4, 3, and 2.

Using standard size 8-1/2- x 11-inch, two-cycle semilog paper, the ordinate or log scale was labeled percent T (1-100) and the linear portion or abscissa was marked off from left to right, one to eight with ten divisions between marks. This axis was labeled NUMBER. Values from the table were plotted and a straight line was drawn connecting the two points corresponding to the values of 6 and 7 and extending to the abscissa axis. After establishment of this line, intermediary points between NUMBERS are taken so that the calibration curve can be drawn more accurately. This was accomplished by choosing an abscissa value between 7 and 8, for example, 7.5, and obtaining the corresponding percent T value from the straight line drawn on the emulsion calibration curve. This value was set equal to a DARK percent T point on the PC and referring back to that curve a corresponding LIGHT percent T value was determined. In the previously established table, place the value 6.5 under NUMBER (6.5 chosen because it is one less than 7.5, the starting value) and the newly obtained LIGHT percent T value under percent T. Continuing, let this new LIGHT percent T value equal that of a DARK percent T point on the PC and obtain a new LIGHT percent T number from the PC. In the table, place 5.5 under NUMBER and the newly obtained LIGHT percent T value under percent T. This procedure is continued until values have been added to the table for NUMBERS 4.5, 3.5, and 2.5. Additional starting

points such as 7.8, 7.6, etc., were chosen and the above example was followed to obtain numerous data points in order to precisely draw the emulsion calibration curve. Typical data points used to establish the preliminary and emulsion calibration curves are shown in Tables 1 and 2.

The emulsion calibration curve just completed was expanded to 15" x 10" semilog paper. The ordinate remained the same scale; however, the abscissa was expanded in such manner that two curves were plotted, one over the other. The NUMBER or abscissa axis scale is enlarged so that the NUMBERS 5 to 8 on the small calibration curve were expanded to full scale for the new curve. Starting with 5, each major division is in tenths and each minor division is equal to 0.02. The top abscissa of the new emulsion calibration curve expands the scale of the small calibration curve region NUMBERS 2 to 5 to full scale. The enlarged curve is completed using values from the table used to construct the small calibration curve.

The last stage of the calibration process involves construction of a rule. A cardboard approximately 28" x 1.5" is used as the starting material and a line is drawn dividing the cardboard in half in the long direction. Abscissa axis NUMBERS 2.0, 3.0, 4.0, 5.0, 6.0, and 7.0 are superimposed on the rule from the expanded calibration curve and marked off from the underside of the long line of the rule. This is done in such manner that approximately abscissa axis NUMBER 7.5 would be the last point at the right hand end of the rule.

With the rule aligned so that NUMBERS 7, 6, and 5 on the graph paper are in line with their corresponding NUMBERS on the rule, a right triangle is used to help superimpose the percent T ordinate values of the large emulsion calibration curve onto the rule and to mark off percent T readings on the upperside of the straight line of the rule. For example, using 2 percent T, one finds the intersection of the ordinate point 2% with the curve, aligns the long straight edge of the triangle at this point while the base of the triangle is parallel with the straight line of the rule, an alignment mark is made, and the NUMBER 2.0 is labeled. In this manner, percent transmissions between 2.0 and approximately 60.0 were marked off. The rule was then shifted so that abscissa axis NUMBERS 2.0, 3.0, and 4.0 on the calibration curve were aligned with the corresponding numbers on the rule and the remainder of the percent T scale was superimposed. The rule has the approximate range 92.0 to 2.0 percent transmission.

Table 1. TYPICAL RAW DARK AND LIGHT DATA POINTS USED TO ESTABLISH THE PRELIMINARY CURVE EQUATION

Point	Dark	Light
1	8.7	24.5
7	5.1	14.7
14	42.1	77.1
15	65.1	91.9
19	11.1	29.1
20	59.5	85.8
30	49.7	81.1
34	30.5	58.1
53	38.1	73.3
73	2.9	6.2

Table 2. TYPICAL LOG PERCENT TRANSMISSION AND NUMBER VALUES FOR ESTABLISHING THE EMULSION CALIBRATION EQUATION

Point	Log Percent Transmission	Percent Transmission
7	0.4047	2.54
6.6	0.5671	3.69
6.2	0.7295	5.36
5.8	0.8917	7.79
5.4	1.0528	11.29
5.2	1.1328	13.58
4.8	1.2907	19.53
4.3	1.4806	30.24
3.7	1.6859	48.52
3.3	1.7974	62.72
2.5	1.9318	85.47

After establishing the rule, one is then ready to interpret densitometer readings of lines on photographic plates in the calibrated wavelength region, corresponding to elements under study. A percent T for the known and for the internal standard is read from the glass emulsion. Going to the rule, and using an ordinary ruler graduated in 24 units to the inch which corresponds to the linear axis spacing on the semilogarithmic paper recommended, the relative intensity difference can be obtained. By setting the zero on the left end of the ruler on the higher percent T, one then notes the point on the ruler which corresponds to the other percent T. That number is the relative intensity for that pair of percent T values. Continuing this process for pairs of known and internal standard lines, an analytical curve can be established on semilog paper with the log axis labeled CONCENTRATION and the linear axis RELATIVE INTENSITY. From this curve unknown concentrations with known relative intensities can be determined. Typical preliminary and emulsion calibration curves and a rule are shown in Figures 1 to 3.

HEWLETT-PACKARD CALCULATOR METHOD

The equation of the preliminary curve can be determined and the curve generated by a plotter by the use of the polynomial regression procedure outlined in the MATH PAC book furnished with the H.P. computer.¹¹ Following this method one inputs a set of data points (Figure 4), and coefficients of a polynomial up to the 20th degree using a least-square fit may be produced. The DARK percent T and LIGHT percent T data used to draw the preliminary curve (see method done by hand) was entered into the H.P. program. Various degrees of the polynomial were




Figure 1. Preliminary curve, exemplifying method of successive determination of DARK $\%T$ and LIGHT $\%T$ values.

11. 9830 A MATH PAC, Model 30, Hewlett-Packard Calculator, v. 1.




Figure 2. Emulsion calibration rule.



Figure 3. Condensed rule.

examined and it was found that a 6th degree gave a sufficiently high correlation coefficient for use. The following equation was used:


$$T_1 = A_0 + A_1 (T_d) + \dots A_6 (T_d)^6$$

where T_1 = percent transmission of light line

T_d = percent transmission of dark line

A_0, A_1, \dots coefficients of polynomial equation.

H. P. CALCULATOR OPERATIONS CHART



T_l = Percent transmission of light line

T_d = Percent transmission of dark line

Figure 4. H.P. calculator operations chart.

The data outputted contained DARK percent T values beginning with the number 1 and incremented in units of 1 to 80 along with corresponding LIGHT percent T values which are the calculated value for each DARK percent T point. Plotting the data, a smooth curve was obtained. This curve was superimposed on a hand-drawn curve from the original data and was found to be essentially identical.

Using the coefficients of the 6th degree polynomial, one generates a table labeled LIGHT percent T and NUMBER which is used to formulate the emulsion calibration curve. Referring to the above equation, let T_d equal 1, solve the above equation for T_l and store this value in the table along with the NUMBER 7. The computer sets $T_d = T_l$ and resolves the equation. In this manner, values in the table for T_l are established for NUMBERS 7, 6, 5, 4, 3, and 2.

A linear curve is assumed between the data points NUMBER = 7, percent T = value found, and NUMBER = 6, percent T = value found. The straight line region is extended to percent T = 1, NUMBER = calculated value from the equation of the straight line. Solving this straight line equation, the computer is then able to take NUMBER values along this line and obtain corresponding T_1 values. For example, using the NUMBER value 7.5, a corresponding T_1 value is obtained from the equation of the straight line. This is set equal to T_d and placed in the previous polynomial equation:

$$T_1 = A_0 + A_1 (T_d) + \dots A_6 (T_d)^6$$

and a corresponding T_1 value is received. In the table T_1 and NUMBER, the T_1 value obtained and the NUMBER 6.5 (one less than 7.5), are stored. In similar fashion, NUMBER values from 6.9 to 2.1 in increments of one tenth are obtained.

The data from the table LIGHT percent T versus NUMBER is then used to generate the coefficients for the equation:

$$\text{NUMBER} = B_0 + B_1 (\log T_1) + \dots B_6 ((\log T_1))^6.$$

A 6th degree polynomial furnished a high correlation coefficient.

Using the above equation, percent T for a known element (obtained from the density of the line on the glass plate emulsion) is substituted for the T_1 value. Solving the equation a value for NUMBER is obtained. The procedure is repeated for a percent T for an internal standard and another NUMBER is obtained. The computer subtracts these two NUMBERS and their difference is called RELATIVE INTENSITY. As outlined in the method done by hand, this process is continued for pairs of known and internal standard lines. An analytical curve is established on semilog paper with the log axis labeled CONCENTRATION and the linear axis RELATIVE INTENSITY. From this curve unknown concentrations with known relative intensities can be determined. Program used is presented in Appendix B.

UNIVAC 1108 COMPUTER METHOD

Appendix C lists the UNIVAC 1108 program. The only input necessary for the program is the initial DARK/LIGHT data read from the emulsion-covered glass plate used for the calibration. Using these values, the computer can calculate coefficients of a polynomial equation to any practical degree required, and from the values obtained, the lowest degree polynomial with good correlation coefficient was chosen. The program then outputted least-square data of DARK/LIGHT points, that is DARK, LIGHT, and LIGHT FITTED, from which the preliminary curve may be drawn. If desired, a plot of this curve can be obtained.

Following the procedure outlined in the Hewlett-Packard Calculator Method, a table labeled LIGHT percent T and NUMBER is generated to formulate the coefficients of the emulsion calibration equation used as the equation of the emulsion calibration curve.

$$\text{NUMBER} = B_0 + B_1 (\log T_1) + \dots B_n (\log T_1)^n$$

where T_1 = percent transmission of light line
 $B_0 \dots B_n$ = coefficients of polynomial equation.

After determining the coefficients of the above equation with the UNIVAC 1108, the equation can be used with many small desk top calculators that do not have the memory needed to solve a 6th degree polynomial. By inserting the log of the percent transmission of an elemental line, one gets the value for NUMBER.

Assuming a calculator is not available, using the above equation, a chart of percent transmission and NUMBER may be generated. Percent T was varied in increments of 0.05 from 2 to 20 and in units of 0.1 from 20 to 90. The left vertical column of the chart lists changes in percent transmission by one unit, for example, from 20 to 21. The top row of the chart varies percent T by either 0.05 or 0.1 depending on the percent transmission range. The body of the chart contains NUMBER values. For example, if an intensity of an element of a standard sample line was measured at 50.5% T and that of an element internal standard line at 20.6% T, one would find the corresponding NUMBER value by proceeding horizontally along the chart from the extreme left column containing the number 50 as the first integer until one reaches the column headed by 0.5. This NUMBER value is then recorded. The procedure is repeated for the other transmission value of 20.6. The difference in NUMBER values is the relative intensity.

Various charts were attempted, including one in which two percent transmissions (for example, from the unknown element line and that of the internal standard) could be found in the first column and first row of the chart. The point of intersection, going horizontally and vertically along the chart from the two percent T values, is the relative intensity and the need of taking the difference in NUMBER for each percent T pair is eliminated. To accomplish this, however, the chart became very voluminous and it was found that the previously described chart was easier to use.

CONCLUSION

A computer program that follows closely the ASTM graphic method has been developed for use on a Hewlett-Packard 9830-A and a UNIVAC 1108 time-sharing computer. For laboratories that do not have continuous access to these computers, one can obtain the coefficients to the equation relating percent transmission to relative intensity and then use the equation on small memory desk top calculators to find relative intensities.

The program used in the UNIVAC 1108 will output a table of percent transmission versus NUMBER which is used to obtain the relative intensity for a pair of percent transmission values.

ACKNOWLEDGMENTS

The author is indebted to Dr. Gary L. Hagnauer for his assistance on the F.P. computer program, to Mrs. Anna M. Hansen for the help on the UNIVAC program, and to Mr. Frederick P. Valente who suggested this timesaving approach.

<p>AD <u>UNCLASSIFIED</u> UNLIMITED DISTRIBUTION</p> <p>Key Words</p> <p>Computer programs Calibrating Spectrographic analysis</p>	<p>AD <u>UNCLASSIFIED</u> UNLIMITED DISTRIBUTION</p> <p>Key Words</p> <p>Computer programs Calibrating Spectrographic analysis</p>	<p>Army Materials and Mechanics Research Center, Watertown, Massachusetts 02172 A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC ANALYSIS - Bernard H. Strauss</p> <p>Technical Report AMMRC TR 77-1, January 1977, 20 pp - illus-tables, D/A Project A1-6-P6350-AM-AM, AMCMS Code 5397.OM.6350</p> <p>A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.</p>	<p>Army Materials and Mechanics Research Center, Watertown, Massachusetts 02172 A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC ANALYSIS - Bernard H. Strauss</p> <p>Technical Report AMMRC TR 77-1, January 1977, 20 pp - illus-tables, D/A Project A1-6-P6350-AM-AM, AMCMS Code 5397.OM.6350</p> <p>A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.</p>
<p>AD <u>UNCLASSIFIED</u> UNLIMITED DISTRIBUTION</p> <p>Key Words</p> <p>Computer programs Calibrating Spectrographic analysis</p>	<p>AD <u>UNCLASSIFIED</u> UNLIMITED DISTRIBUTION</p> <p>Key Words</p> <p>Computer programs Calibrating Spectrographic analysis</p>	<p>Army Materials and Mechanics Research Center, Watertown, Massachusetts 02172 A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC ANALYSIS - Bernard H. Strauss</p> <p>Technical Report AMMRC TR 77-1, January 1977, 20 pp - illus-tables, D/A Project A1-6-P6350-AM-AM, AMCMS Code 5397.OM.6350</p> <p>A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.</p>	<p>Army Materials and Mechanics Research Center, Watertown, Massachusetts 02172 A COMPUTER METHOD FOR THE CALIBRATION OF A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC ANALYSIS - Bernard H. Strauss</p> <p>Technical Report AMMRC TR 77-1, January 1977, 20 pp - illus-tables, D/A Project A1-6-P6350-AM-AM, AMCMS Code 5397.OM.6350</p> <p>A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.</p>

Army Materials and Mechanics Research Center,
Watertown, Massachusetts 02172
A COMPUTER METHOD FOR THE CALIBRATION OF
A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC
ANALYSIS - Bernard H. Strauss

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION
Key Words

Technical Report AMMRC TR 77-1, January 1977, 20 pp -
illus-tables, D/A Project A1-6-P6350-AW-AW,
AMCMS Code 5397.OM.6350

A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.

Army Materials and Mechanics Research Center,
Watertown, Massachusetts 02172
A COMPUTER METHOD FOR THE CALIBRATION OF
A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC
ANALYSIS - Bernard H. Strauss

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION
Key Words

Technical Report AMMRC TR 77-1, January 1977, 20 pp -
illus-tables, D/A Project A1-6-P6350-AW-AW,
AMCMS Code 5397.OM.6350

A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.

Army Materials and Mechanics Research Center,
Watertown, Massachusetts 02172
A COMPUTER METHOD FOR THE CALIBRATION OF
A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC
ANALYSIS - Bernard H. Strauss

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION
Key Words

Technical Report AMMRC TR 77-1, January 1977, 20 pp -
illus-tables, D/A Project A1-6-P6350-AW-AW,
AMCMS Code 5397.OM.6350

A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.

Army Materials and Mechanics Research Center,
Watertown, Massachusetts 02172
A COMPUTER METHOD FOR THE CALIBRATION OF
A PHOTOGRAPHIC EMULSION FOR SPECTROGRAPHIC
ANALYSIS - Bernard H. Strauss

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION
Key Words

Technical Report AMMRC TR 77-1, January 1977, 20 pp -
illus-tables, D/A Project A1-6-P6350-AW-AW,
AMCMS Code 5397.OM.6350

A computer method to calibrate emulsion-covered glass plates used to interpret densities of lines for quantitative spectrographic work is reported. This method follows a slightly modified ASTM graphic method, which is also described, and allows one to simply input dark and light data points to establish preliminary and emulsion calibration curves. After the curves are established, by inputting percent transmission data pairs for a known or unknown concentration line and an internal standard line, one obtains the relative intensity. The computer method is used on a Hewlett-Packard Calculator 9830-A and a UNIVAC 1108 time-sharing computer. Programs for both computers are included in this report.

APPENDIX A. SPECTROGRAPHIC PARAMETERS

Equipment

Baird 3-meter spectrograph, 2.75 Å/mm reciprocal linear dispersion
Jarrell-Ash Varisource
Jarrell-Ash Microphotometer

Instrumental Parameters

Spark-18,000 V. A. C., 3 breaks per half-cycle, 6 to 6.5 amperes
Analytical gap-3 mm
Auxiliary gap-6 mm
Sample electrode-Ultra Carbon 106
Counter electrode-Ultra Carbon 105 U
Pre-spark duration-20 sec.

Photographic Parameters

Plates-Eastman Kodak SA 1
Developer-Eastman Kodak D19
Developing time- 3.5 min.
Developing temperature 20°C

Spectral range of calibration 2750-3450 Å
Spectral lines read in range of 3000-3200 Å

APPENDIX B. HEWLETT-PACKARD CALCULATOR PROGRAM FOR GLASS PLATE EMULSION CALIBRATION

f₀

```

10  For I=1 to 22
20  C [I] = B [I] = 0
30  Next I
40  For I = 23 to 253
50  C [I] = 0
60  Next I
70  B [I] = I
80  W = N = S1 = S2 = S3 = S4 = S5 = I1 = 0
90  Disp " Max Degree Polynomial ";
100 Input D2
110 Print " Max Degree Polynomial = " D2
120 D1 = D2
130 Print
140 End

```

f₁

```

10  Format F 4.0, 2 F 12.4
20  Print "Enter Prelim Calibrn Data"
30  If N ≠ 0 then 70
40  Print "No." Tab 10" Dark "Tab 22 "Light"
50  I1 = 0
60  I1 = I1 + 1
70  Disp "Dark, Light =";
80  Input G[I1], H [I1]
90  Y = H [I1]
100 B [2] = G[I1]
110 Write (15, 10) N+1, G [I1], H [I1]
120 Y = FN X 1
130 Go to 60
140 End

```

f₂

```

10  If W = 0 then 40
20  Disp " Not allowed"

```

```

30  End
40  Disp " Wrong No., Dark, Light =";
50  Input I, Y, B [2]
60  G [I] = Y
70  H [I] = B [2]
80  Print "Delete: Dark = "G[I]" Light = "H [I]
90  Y = FNX (-1)
100 Disp
110 End

```

f₃

```

10  I1 = 0
20  Print " Polynomial Degree = "D1
30  Print "Relative I Log T Light "
40  X2 = 1
50  X = 1
60  X1 = 8
70  If X1 <= 2 then 210
80  Y = B[D1 +1]
90  For J = D1 to 1 Step -1
100 Y = Y * X + B[J]
110 Next J
120 X=Y
130 If X>100 then 210
140 X1 = X1-1
150 If X1 <= 2 then 210
160 I1 = I1 + 1
170 H [I1] = X1
180 G [I1] = L G T X
190 Print X1, G [I1], X
200 Go to 70
210 X2 = X2 - 0.1
220 If X2 = 0 then 290
230 I1 = I1 + 1

```

```

240  G [11] = (G [2] - G [3]) * X2 + G [2]
250  X1 = 6 + X2
260  X = 10*G [11]
270  H [11] = X1
280  Go to 190
290  Y = FN Y1
300  End

f4
10  Def FNX (Z)
20  For I = 2 to D2
30  B [I+1] = B [I] * B [2]
40  Next I
50  B [D2 + 2] = Y
60  R = 0
70  For I = 1 to D2 + 2
80  For J = 1 to D2 + 2
90  R = R+1
100 C [R] = C [R] + B [I] * B[J] * Z
110 Next J
120 Next I
130 S1 = S1 + B [2] * Z
140 S2 = S2 + B [2]*2 * Z
150 S3 = S3 + Y * Z
160 S4 = S4 + Y *Y * Z
170 S5 = S5 + B [2] * Y * Z
180 N = N + Z
190 Return 0

f5
10  Disp " New Polynomial Degree =";
20  Input D1
30  Print " Degree of Polynomial =" D1
40  End

f6
10  If N <= D2 - W then 340
20  If D1 <= D2 - W then 50
30  Disp " Max Deg ="; D2 - W
40  End
50  If W = 0 then 330
60  T = 0
70  For I = 1 to D1 + 1
80  B [I] = 0
90  For J = 1 to D1 - I + 2
100 R = ( I + J - 1) * D2 + 2 - 0.5 * ( I + J )
110 B [I] = B [I] + C [ T + J] * C [R]
120 Next J
130 T = I * ( D2 + (3-1) / 2)
140 Next I
150 R1 = 0
160 For I = 2 to D1 + 1
170 R1 = R1 + C [I*(D2 + (3-1) / 2)+2]
180 Next I
190 T0 = C [(D2 + 1) * (D2+2) / 2]
200 T0 = T0 - C [D2 + 1]+2
210 Print
220 Print " Coefficients"
230 Print
240 Format F3.0, E 12.4
250 For I = 1 to D1 + 1
260 Write (15, 240) " B ("I-1") = "B[I]
270 Next I
280 Print
290 Print
300 Print "R Square = " R1/T0
310 Print
320 End
330 If N > D2 then 360
340 Disp " Not Enough Points"
350 End
360 P = W = 1

```



```

370 D2 = D2 + 1
380 For J = 1 to D2
390 Print C [P]
400 C [P] = SQR C [P]
410 For I = 1 to D2 - J+1
420 C [P+I] = C [P+I] / C [P]
430 Next I
440 R = P + I
450 S = R
460 For L = 1 TO D2 - J
470 P = P + 1
480 For M = 1 to D2 + 2 - J - L
490 C [ R+M - 1] + C [ R+M-1] - C [P] * C [P+M-1] f8
500 Next M
510 R = R + M - 1
520 Next L
530 P = S
540 Next J
550 T = ( D2 + 1) * ( D2 + 2) / 2
560 For I = 1 to D2 - 1
570 T = T - I - I
580 C [T] = 1 / C [T]
590 For J = 1 to D2 - 1
600 P = D2 + 1 - I - J
610 P = P * ( D2 + 1 - (P-1) / 2) - I
620 R = P-J
630 S = 0
640 U = I + J + 1
650 Y = P
660 For K = 1 to J
670 Y = Y + U - K
680 S = S - C[R+K] * C [Y]
690 Next K
700 C [P] = S / C [R]

710 Next J
720 Next I
730 C [1] = 1/ C [1]
740 Go to 60

f7
10 B [22] = D1
20 Print " Store X = F ( Log T) Calibration "
30 Disp " Store File # = " ;
40 Input I
50 Print " Store File # " I
60 Store Data I, B
70 End

10 Print " Diff X Evaluation"
20 Disp " T (STD), T = ";
30 Input T1, T2
40 Print " T (STD) = "T1
50 Print " T = " T2
60 T1 = LGTT1
70 T2 = LGTT2
80 X1 = X2 = B [D1 + i]
90 For J = D1 to 1 Step - 1
100 X1 = X1 * T1 + B [J]
110 X2 = X2 * T2 + B [J]
120 Next J
130 Print " X (STD) = " x1
140 Print " X = " X2
150 Print " X (STD) - X = "X1-X2
160 Go to 20
170 End

f9
10 Disp "Re-plot Axis";
20 Input Z9
30 If Z9 = 0 then 390

```

```

40    Disp "X (Max) = ";
50    Input B2
60    Disp "X (Min) = ";
70    Input A2
80    Disp "Y(Max) = ";
90    Input B1
100   Disp "Y (Min) = " ;
110   Input A1
120   C2 = (B2 - A2) / 10
130   C1 = (B1 - A1) / 10
140   Scale A2, B2 + 2*C2, A1, B1, + 2*C1
150   Plot A2 + (B2 - A2) / 2.2, A1 + C1 / 4
160   Label (*, 3, 2, 0, 2/3)
170   Letter
180   Plot A2 + C2/4, A1 + (B1 - A1)/2.2
190   Label (*, 3, 2, P1/ 2, 2/3)
200   Letter
220   Offset C2, C1
230   XAXIS A1, C2, A2, B2
240   YAXIS B2, C1, A1, B1
250   YAXIS A2, C1, A1, B1
260   XAXIS B1, C2, A2, B2
270   Label (*, 2, 1.7, 0, 2/3)
280   For Y = A1 to B1 Step 2*C1
290   Plot A2, Y, 1
300   CPL0T -5.4, -0.34
310   Label (380, 2, 1.7, 0, 2/3) Y
320   Next Y
330   For X = A2 to B2 step 2*C2
340   Plot X, A1, 1
350   CPL0T -3.8, -1
360   Label (380, 2, 1.7, 0, 2/3) X
370   Next X
380   Format F 5.1

```

```

390   For X = A2 to B2 Step C2/10
400   Y = B [D1 + 1]
410   For J = D1 to 1 Step -1
420   Y = Y * X + B [J]
430   Next J
440   Plot X,Y
450   Next X
460   Disp "Plot Data Points";
470   Input Z9
480   If Z9 = 0 then 550
490   PEN
500   For K = 1 to 11
510   Plot G[K], H [K]
520   CPL0T -0.3, -0.3
530   Label (*) "0"
540   Next K
550   PEN
560   Plot -1.2 * C2, -C1
570   PEN
580   Go to 10

f11
10   Format 2F 12.4
20   Print Tab 9 "X Dark "Tab 21" Y Light"
30   For I = 1 to 79 Step 2
40   Y = B [D1 + 1]
50   For J = D1 to 1 Step -1
60   Y = Y* I + B [J]
70   Next J
80   Write (15, 10) I, Y
90   Next I
100  For I = 80 to 100 Step 1
110  Y = B [Di + 1]
120  For J = D1 to 1 Step -1
130  Y = Y * I + B[J]

```



```

140 Next J
150 Write (15, 10) I, Y
160 Next I
170 End

f13

10 Format 2 F 12.4
20 Print Tab 9 "X" Tab 21 "T"
30 For I = 0.3 to 2 Step 0.05
40 Y = B [D1 + I]
50 For J = D1 to 1 Step -1
60 Y = Y * I + B [J]
70 Next J
80 Write (15, 10) Y, 10+I
90 Next I
100 End

f14

10 Def FNY(Z)
20 Disp "Max Degree Polynomial";
30 Input D2
40 Print "Degree Polynomial = "D2
50 D1 = D2
60 Print I1
70 For I = 1 to 22
80 C [I] = B[I] = 0
90 Next I
100 For I = 23 to 253
110 C [I] = 0
120 Next I
130 I3 = 0
140 B [1] = 1
150 W = N = S1 = S2 = S3 = S4 = S5 = 0
160 I3 = I3 + 1
170 B [2] = G [I3]
180 Y = H [I3]

```

```

190 Y = FN X1
200 Print I3, G[I3], H [I3]
210 IF I3 = 11 - 1 then 230
220 Go to 160
230 End

f17

10 Print "Load X = F (LOG T) Calibration"
20 Disp "Load File # = ";
30 Input I
40 Print "Load File # "I
50 Load Data I, B
60 D1 = D2 = B [22]
70 Print "Degree Polynomial = "D1
80 End

f18

10 Disp "Store Data, File # = ";
20 Input I
30 Store Data I
40 End

f19

10 Disp "Load Data, File # = ";
20 Input I
30 Load Data I
40 Format F 4.0, 2 F 12.4
50 Print "No." Tab 10 "Dark" Tab 22 "Light"
60 I1 = 0
70 I1 = I1 + 1
80 Y = H [I1]
90 B[2] = G [I1]
100 Write (15, 40) N + 1, G [I1], H[I1]
110 Y = FN X1
120 Go to 70
130 END

```

APPENDIX C. UNIVAC 1108 PROGRAM

```

        DIMENSION X(1200), Y(1200), W(200), T1(200), T2(200), T3(200),
1      C(11), ALPHA(10), BETA(10),
        DIMENSION S(4), YF(400), A(14)
105     FORMAT (/// 29H OUTPUT FROM SUBROUTINE FITY.//7H IND 3 = I 3///3X
        2H 1X, 18X 2HY, 14X 9H Y FITTED // (1H 1P3E 20.7))
103     FORMAT (///31H OUTPUT FROM SUBROUTINE ORTHLS.//7H IND 1 = 13//2X
1      2HI, 13X 2H C, 16X 6H ALPHA, 14X 5H BETA/1H I4, 1PE 20.7/
2      (1H I4, 1P3E 20.7))
104     FORMAT (/// 30H OUTPUT FROM SUBROUTINE COEFS.//7H IND 2 = 13//13X
1      2HI, 17X 5H A(I)//(1H I14, 1PE 26.7))
        P = 3.141592654
        READ 1,N,K
        READ 1, NUM, NUMA, DUM
1      FORMAT (I5, I5, F10.0)
        DO 10 I = 1, N
        READ 2, X(I), Y(I)
10     CONTINUE
2      FORMAT (F10.0, F10.0)
101    L = 0
        J = 0
        CALL ORTHLS (X, Y, W, N, L, J, C, ALPHA, BETA, K, T1, T2, T3, IND1)
        IT = 0
        WRITE (6, 103) IND 1, II, C(1), (II, C(II + 1), ALPHA (II), BETA (II),
1      II = 1, K)
        DO 15 MM = 1,K
        KF = MM
        CALL FITY (X, N, J, C, ALPHA, BETA, KF, YF, T1, T2, IND 3)
        WRITE (6, 105) IND 3, (X(I), Y(I), YF(I), I = 1, N)
        CALL COEFS (J, C, ALPHA, BETA, KF, A, T1, T3, IND 2)
        IT = 0
        WRITE (6, 104) IND2, II, A(1), II, A(II + 1),
        II = 1, KF)
82     FORMAT (1H1)
81     FORMAT (1H 3E 15.6)
        MB = MM +1
15     CONTINUE
        PRINT 82
        DO 80 IA = 1, N
        ANS = A(1)
        DO 75 II = 2, MB
        FIB = II -1
75     ANS = ANS + A(II) *X(IA)* FIB
        IF (IT. EQ.3) X(IA) = EXP (X(IA))
80     PRINT 81, X(IA), Y(IA), ANS

```

```

IF, (IT.EQ.3) GO TO 102
PRINT 82
AB = 1.0
SAVE = SOLVE (A, AB, MB)
PRINT 85, SAVE, AB
SAVE = SOLVE (A, SAVE, MB)
PRINT 85, SAVE 1, SAVE
CON = ALOG (SAVE)
CON 1 = ALOG (SAVE 1)
D = CON1-CON
Y(1) = ALOG (SAVE)
X(1) = 7
S(1) = SAVE
DD = D/10
DO 650 I = 2, 10
S(I) = EXP (CON + DD * (I - 1))
X(I) = 7. -(I-1)*.1
650 Y(I) = ALOG (S(I))
L = 10
DO 651 J = 1, 10
AB = S(J)
DO 652 JJ = 1,4
L = L + 1
ANS = SOLVE (A, AB, MB)
Y(L) = ALOG (ANS)
X(L) = X(J) -1.0 * JJ
IF (ANS. GE. 100.) GO TO 653
PRINT 85, ANS, AB, X(L)
652 AB = ANS
GO TO 651
653 L = L-1
651 CONTINUE
IT = 3
N = L
DO 107 J = 1, N
TEM = X(J)
X(J) = Y(J)
107 Y(J) = TEM
GO TO 101
102 CONTINUE
NUMB = NUM + 10
NUMC = NUMA + 1
DO 111 I = 1, NUMB
X(I) = DUM
IF (I. GT. NUMC) X(I) = X(1) + .05

```

```

AB = ALOG (X(I))
Y(I) = SOLVE (A, AB, MB)
111 DUM = X(1) + .05
PRINT 17
WRITE (6, 18)
DO 40 I = 1, NUMA, 20
JJ = I + 19
40 WRITE (6, 19) X(I), (Y(J), J = I, JJ)
NUMB = JJ + 1
IPRINT = NUMB + 500
PRINT 17
WRITE (6, 13)
DO 150 I = NUMB, NUM, 10
JJ = I + 9
IF(I. GT. IPRINT) GO TO 16
GO TO 150
16 IPRINT = I + 500
PRINT 17
PRINT 13
150 WRITE (6,11) X(I), (Y(J), J = I, JJ)
19 FORMAT (1H 21F6.2)
18 FORMAT (16 X, '.05' 3X, '.10' 3X, '.15' 3X, '.20' 3X, '.25'
3X, '.30',
1 3X, '.35' 3X, '.40' 3X, '.45' 3X, '.50' 3X, '.55' 3X, '.60' 3X, '.65'
2 3X, '.70' 3X, '.75' 3X, '.80' 3X, '.85' 3X, '.90' 3X, '.95'
11 FORMAT (1H F6.2, 10F8.2)
13 FORMAT (1H 19X '.1', 6X '.2', 6X '.3', 6X '.4', 6X '.5', 6X '.6',
6X '.7', 6X '.8', 6X '.9')
17 FORMAT (1H1)
85 FORMAT (1H 'Y 1 = ' F13.8, 5X, 'X1 = ' F13.8, E14.8)
STOP
END
LEGEND
N = No. OF EXPERIMENTAL POINTS
K = MAXIMUM ORDER OF EQUATION FOR POLYNOMIAL LEAST SQUARES FIT.
NUM = TOTAL NUMBER OF PERCENT TRANSMISSION vs RELATIVE INTENSITY VALUES
IN TABLE
NUMA = TOTAL NUMBER OF PERCENT TRANSMISSION VALUES OF INCREMENT .05 IN TABLE
DUM = INITIAL VALUE OF PERCENT TRANSMISSION PRINTED IN TABLE (2.0 USED IN
THIS WORK)
X(I) = EXPERIMENTAL DARK PERCENT TRANSMISSION
Y(I) = EXPERIMENTAL LIGHT PERCENT TRANSMISSION

```


DISTRIBUTION LIST

No. of Copies	To	No. of Copies	To
2	Metals and Ceramics Information Center, Battelle Columbus Laboratories, 505 King Avenue, Columbus, Ohio 43201	1	Commander, Edgewood Arsenal, Edgewood, Maryland 21010
12	Commander, Defense Documentation Center, Cameron Station, Building 5, 5010 Duke Street, Alexandria, Virginia 22314	1	ATTN: SAREA-TSP
	Office of Chief of Research and Development, Department of the Army, Washington, D. C. 20310	1	SAREA-QAE
1	ATTN: CRDPES	1	SAREA-QAI
1	Commander, Army Research Office, P. O. Box 12211, Research Triangle Park, North Carolina 27709	1	SAREA-QAP
	Commander, U. S. Army Material Development and Readiness Command, 5001 Eisenhower Avenue, Alexandria, Virginia 22333	1	SAREA-QAIP
1	ATTN: DRCQA		Commander, U. S. Army Foreign Science and Technology Center, 220 7th Street, N. E., Charlottesville, Virginia 22901
2	DRCQA-E	1	ATTN: DRXST-SD2
2	DRCQA-P		Commander, Frankford Arsenal, Philadelphia, Pennsylvania 19137
1	DRCLDC, Mr. R. Zentner	1	ATTN: SARFA-P3300
1	DRCLRL	1	SARFA-C2500
1	DREDE-EA	1	SARFA-Q1000
1	DRCDE-E	2	SARFA-Q2100
1	DRCRP-OIP	1	SARFA-N3100-202-1, Mr. E. Roffman
	Commander, U. S. Army Electronics Command, Fort Monmouth, New Jersey 07703	1	SARFA-Q6120-64-1, Mr. W. Shebest
2	ATTN: DRSEL-CB	1	SARFA-Q6130-64-1, Mr. G. Norwitz
1	DRSEL-PP-PQ	1	SARFA-A2000
1	DRSEL-WM	2	SARFA-F6000
1	DRSEL-RD-GT	1	SARFA-L300, Mr. J. Corrie
1	DRSEL-PA-C		Commander, Harry Diamond Laboratories, 2800 Powder Mill Road, Adelphi, Maryland 20783
	Commander, U. S. Army Missile Command, Redstone Arsenal, Alabama 35899	2	ATTN: DRXDO-EDE
2	ATTN: DRSMI-RBLD, Redstone Scientific Information Center		Commander, Picatinny Arsenal, Dover, New Jersey 07801
1	DRSMI-RSM, Mr. E. J. Wheelahan	1	ATTN: SARPA-RT-S
1	DRSMI-RTR, Mr. H. T. Lawson	1	SARPA-VA6, Mr. H. DeFazio
2	DRSMI-Q	1	SARPA-VC2, Mr. T. M. Roach, Jr.
1	DRSMI-M	1	SARPA-VG, Mr. A. Clear
	Commander, U. S. Army Mobility Equipment Research and Development Command, 4300 Goodfellow Boulevard, St. Louis, Missouri 63120	1	SARPA-ND 1, Mr. D. Stein
1	ATTN: DRSMQ	1	SARPA-FR-M-D, A. M. Anzalone, Bldg. 176
1	DRSMQ-QP		Commander, Rock Island Arsenal, Rock Island, Illinois 61201
1	DRSMQ-QR	1	ATTN: SARRI-L, Dr. Royce Beckett
1	DRSMQ-QE	1	SARRI-LEQ
1	DRSMQ-M	2	SARRI-LRM
1	DRSMQ-P		Director, U. S. Army Production Equipment Agency, Rock Island Arsenal, Rock Island, Illinois 61201
1	DRSMQ-R	1	ATTN: DRXDE-MT
	Commander, U. S. Army Armament Command, Rock Island, Illinois 61201		Commander, U. S. Army Aeronautical Depot Maintenance Center, Corpus Christi, Texas 78419
1	ATTN: DR SAR-SC, Dr. C. M. Hudson	1	ATTN: SSMAC-Q
1	DR SAR-PPW-PB, Mr. Francis X. Walter		Commander, U. S. Army Ammunition Procurement and Supply Agency, Joliet, Illinois 60436
2	Technical Library	2	ATTN: SMUAT-E
2	Commander, U. S. Army Natick Research and Development Command, Kansas Street, Natick, Massachusetts 01760		Commander, U. S. Army Mobility Equipment Research and Development Command, Fort Belvoir, Virginia 22060
1	ATTN: DRXNM-QE	1	ATTN: Technical Documents Center, Bldg. 315
1	DRXNM-GE	1	DRXFB-P
	Commander, U. S. Army Tank-Automotive Research and Development Command, 28251 Van Dyke Avenue, Warren, Michigan 48090	1	DRXFB-M
1	ATTN: DRDTA-RXMT, Mr. J. Dudzinski	1	DRXFB-X
1	DRDTA-RXMT, Mr. F. Lemmer	1	DRXFB-A
1	DRDTA-RXMT, Mr. W. Wulf	1	DRXFB
1	DRDTA-RGD, Mr. O. Renius	1	DRXFB-B
2	DRDTA-QEE, Mr. P. Duika	1	DRXFB-H
	Commander, U. S. Army Test and Evaluation Command, Aberdeen Proving Ground, Maryland 21005	1	DRXFB-J
2	ATTN: DRSTE-TA-A	1	DRXFB-F
	Commander, Aberdeen Proving Ground, Maryland 21005	1	DRXFB-MM
1	ATTN: STEAP-MT, Mr. J. M. McKinley	1	DRXFB-Q
1	STEAP-TL	1	DRXFB-QQ
		1	DRXFB-QE, Mr. Jacob K. Mauzy
		1	DRXFB-MW, Dr. J. W. Bond
		2	Commander, Watervliet Arsenal, Watervliet, New York 12189
		1	ATTN: SARWV-QA, Mr. J. Miller
		1	SARWV-QA, Quality Assurance Office
			Commander, Anniston Army Depot, Anniston, Alabama 36202
		1	ATTN: AMXAN-QA
			Commander, Atlanta Army Depot, Forest Park, Georgia 30050
		1	ATTN: AMXAT-CSQ

No. of Copies	To
1	Commander, Letterkenny Army Depot, Chambersburg, Pennsylvania 17201
1	ATTN: AMXLE-CQ
1	ATTN: AMXLE-NSQ
1	Commander, Lexington-Bluegrass Army Depot, Lexington, Kentucky 40507
1	ATTN: AMXLX-QA
1	Commander, New Cumberland Army Depot, New Cumberland, Pennsylvania 17070
1	ATTN: AMXNC-256
1	Commander, Pueblo Army Depot, Pueblo, Colorado 81001
1	ATTN: AMXPU-BF
1	Commander, Red River Army Depot, Texarkana, Texas 75502
1	ATTN: AMXRR-QAO
1	Commander, Sacramento Army Depot, Sacramento, California 95801
1	ATTN: AMXSA-QA
1	Commander, Savanna Army Depot, Savanna, Illinois 61074
1	ATTN: AMXSV-QAO
1	Commander, Seneca Army Depot, Romulus, New York 14541
1	ATTN: AMXSE-AXI
1	Commander, Sharpe Army Depot, Lathrop, California 95330
1	ATTN: AMXSH-CQ
1	Commander, Sierra Army Depot, Herlong, California 96113
1	ATTN: AMXSI-QA

No. of Copies	To
1	Commander, Tobyhanna Army Depot, Tobyhanna, Pennsylvania 18466
1	ATTN: AMXTQ-Q
1	Commander, Tooele Army Depot, Tooele, Utah 84074
1	ATTN: AMXTE-QAD
1	Commander, Umatilla Army Depot, Hermiston, Oregon 97838
1	ATTN: AMXUM-QA
1	Chief, Bureau of Naval Weapons, Department of the Navy, Washington, D. C. 20390
1	Chief, Bureau of Ships, Department of the Navy, Washington, D. C. 20315
1	Commander, Wright Air Development Division, Wright-Patterson Air Force Base, Ohio 45433
2	ATTN: ASRC
1	AFML/MBC/Mr. Stanley Schulman
1	Commander, U. S. Army Aviation Systems Command, St. Louis, Missouri 63166
1	ATTN: SRSV-R-R
1	SRSV-R-EGE
1	SRSV-A-L
1	SRSV-LE
1	SRSV-A-LV
1	SRSV-A-V
1	Director, Army Materials and Mechanics Research Center, Watertown, Massachusetts 02172
2	ATTN: DRXMR-PL
1	DRXMR-AG
2	Authors